

Product Datasheet

Bovine CCL4/C-C motif chemokine 4 ELISA Kit (orb1662244)

| | |
|----------------------------|---|
| Description | Bovine CCL4/C-C motif chemokine 4 ELISA Kit |
| Reactivity | Bovine |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 15.6-1000 pg/mL. Sensitivity: 7.8 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: CCL4. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| Expiration Date | Please enquire. |