

Product Datasheet

Human RAR Related Orphan Receptor Alpha (RORa) ELISA Kit (orb1662434)

Description	Human RAR Related Orphan Receptor Alpha (RORa) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.056 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.