




Product Datasheet

Pbk antibody (orb166244)

Description	Rabbit polyclonal antibody to Pbk (N-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	Pbk
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36745
Uniprot ID	Q9JJ78
NCBI	NP_075698.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of A431 cells u...



Immunohistochemical analysis of formalin...