

Product Datasheet

Human JNK1/MAPK8 Associated Membrane Protein (JKAMP) ELISA Kit (orb1662632)

Description	Human JNK1/MAPK8 Associated Membrane Protein (JKAMP) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.781-50 ng/mL. Sensitivity: 0.35 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.