

Product Datasheet

Human Wingless Type MMTV Integration Site Family, Member 11 (WNT11) ELISA Kit (orb1662807)

Description	Human Wingless Type MMTV Integration Site Family, Member 11 (WNT11) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.058 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.