

Product Datasheet

Human Von Willebrand Factor Cleaving Protease (vWFCP) ELISA Kit (orb1662824)

Description Human Von Willebrand Factor Cleaving Protease (vWFCP) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 62.5-4000 pg/mL. Sensitivity: 21.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.

Expiration Date Please enquire.