

Product Datasheet

General Vitamin D2 (VD2) ELISA Kit (orb1662887)

Description	General Vitamin D2 (VD2) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 4.68-300 ng/mL. Sensitivity: 2.2 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.