

## Product Datasheet

# Human Ubiquitin Cross Reactive Protein (UCRP) ELISA Kit (orb1662939)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human Ubiquitin Cross Reactive Protein (UCRP) ELISA Kit   |
| <b>Reactivity</b>          | Human   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detect Range: 31.2-2000 pg/mL. Sensitivity: 12.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.   |