

Product Datasheet

Mouse Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) ELISA Kit (orb1662958)

Description	Mouse Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 25 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.