

Product Datasheet

Human Ubiquitin Conjugating Enzyme E2C (UBE2C) ELISA Kit (orb1662973)

Description	Human Ubiquitin Conjugating Enzyme E2C (UBE2C) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.053 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.