

Product Datasheet

Mouse Thyrotropin Releasing Hormone (TRH) ELISA Kit (orb1663017)

Description	Mouse Thyrotropin Releasing Hormone (TRH) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.81-500 pg/mL. Sensitivity: 2.77 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.