

Product Datasheet

Porcine Plasminogen Activator, Tissue (tPA) ELISA Kit (orb1663066)

Description	Porcine Plasminogen Activator, Tissue (tPA) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.781-50 ng/mL. Sensitivity: 0.32 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: plasma or other biological fluids.
Expiration Date	Please enquire.