

Product Datasheet

Canine Tumor Necrosis Factor Beta (TNFb) ELISA Kit (orb1663145)

Description	Canine Tumor Necrosis Factor Beta (TNFb) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 ng/mL. Sensitivity: 3.6 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.