

Product Datasheet

Bovine Tumor Necrosis Factor Beta (TNFb) ELISA Kit (orb1663146)

Description	Bovine Tumor Necrosis Factor Beta (TNFb) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.041 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.