

## Product Datasheet

# Porcine Thrombomodulin (TM) ELISA Kit (orb1663169)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Porcine Thrombomodulin (TM) ELISA Kit   |
| <b>Reactivity</b>          | Porcine   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detect Range: 0.781-50 ng/mL.<br>Sensitivity: 0.29 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.<br>Sample Volume: 100 µL. Sample Type: plasma. |
| <b>Expiration Date</b>     | Please enquire.   |