

Product Datasheet

Canine Thrombomodulin (TM) ELISA Kit (orb1663175)

Description	Canine Thrombomodulin (TM) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 5.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: plasma.
Expiration Date	Please enquire.