

## Product Datasheet

# Mouse Toll Like Receptor 6 (TLR6) ELISA Kit (orb1663180)

<b>Description</b>	Mouse Toll Like Receptor 6 (TLR6) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.125 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Expiration Date</b>	Please enquire.