

Product Datasheet

Human Toll Like Receptor 2 (TLR2) ELISA Kit (orb1663193)

Description	Human Toll Like Receptor 2 (TLR2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.113 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, saliva, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.