

Product Datasheet

Porcine Tissue Inhibitors Of Metalloproteinase 1 (TIMP1) ELISA Kit (orb1663224)

Description	Porcine Tissue Inhibitors Of Metalloproteinase 1 (TIMP1) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 13.6 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.