

## Product Datasheet

# Porcine Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit (orb1663290)

<b>Description</b>	Porcine Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 31.2-2000 pg/mL. Sensitivity: 10.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, platelet-poor plasma, tissue homogenates, cell culture supernates or other biological fluids.
<b>Expiration Date</b>	Please enquire.