

Product Datasheet

Rat Transferrin Receptor (TFR) ELISA Kit (orb1663304)

Description Rat Transferrin Receptor (TFR) ELISA Kit

Reactivity Rat

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.625-40 ng/mL. Sensitivity: 0.236 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.

Expiration Date Please enquire.