

Product Datasheet

Human Thymosin Beta 4 (Tb4) ELISA Kit (orb1663350)

Description Human Thymosin Beta 4 (Tb4) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 3.12-200 ng/mL. Sensitivity: 1.25 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 μ L. Sample Type: serum, plasma or other biological fluids.

Expiration Date Please enquire.