

Product Datasheet

Human Suppressor Of Zeste 12 Homolog (SUZ12) ELISA Kit (orb1663379)

Description	Human Suppressor Of Zeste 12 Homolog (SUZ12) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.113 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.