

Product Datasheet

Human Steroid Sulfatase Isozyme S (STS) ELISA Kit (orb1663396)

Description	Human Steroid Sulfatase Isozyme S (STS) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 4.68-300 ng/mL. Sensitivity: 1.9 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.