

Product Datasheet

Human Surfactant Associated Protein B (SPB) ELISA Kit (orb1663458)

Description	Human Surfactant Associated Protein B (SPB) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.54 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.