

Product Datasheet

CB1 antibody (orb166347)

Description	Rabbit polyclonal antibody to CB1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CNR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52858
Uniprot ID	P21554
NCBI	NP_001153698.1 , NP_001153730.1 , NP_001153731.1 , NP_149421.2 , NP_057167.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

