

Product Datasheet

Human Solute Carrier Family 39, Member 6 (SLC39A6) ELISA Kit (orb1663524)

Description	Human Solute Carrier Family 39, Member 6 (SLC39A6) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.114 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.