

Product Datasheet

Human Signaling Lymphocytic Activation Molecule Family, Member 7 (SLAMF7) ELISA Kit (orb1663540)

Description	Human Signaling Lymphocytic Activation Molecule Family, Member 7 (SLAMF7) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.053 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.