

Product Datasheet

Human Sialic Acid Binding Ig Like Lectin 10 (SIGLEC10) ELISA Kit (orb1663564)

Description	Human Sialic Acid Binding Ig Like Lectin 10 (SIGLEC10) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.062 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.