

Product Datasheet

Human Sialic Acid Acetylerase (SIAE) ELISA Kit (orb1663567)

Description	Human Sialic Acid Acetylerase (SIAE) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.054 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.