

Product Datasheet

Canine Selectin, Endothelium (SELE) ELISA Kit (orb1663640)

Description	Canine Selectin, Endothelium (SELE) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 ng/mL. Sensitivity: 2.38 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.