

Product Datasheet

General Salmon Calcitonin (SCT) ELISA Kit (orb1663656)

Description	General Salmon Calcitonin (SCT) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.085 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 μ L. Sample Type: biological agents.
Expiration Date	Please enquire.