

Product Datasheet

General Reverse Triiodothyronine (rT3) ELISA Kit (orb1663727)

| | |
|----------------------------|---|
| Description | General Reverse Triiodothyronine (rT3) ELISA Kit |
| Reactivity | All |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 93.7-6000 pg/mL. Sensitivity: 24.32 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids. |
| Expiration Date | Please enquire. |