

## Product Datasheet

### Porcine Relaxin (RLN) ELISA Kit (orb1663769)

<b>Description</b>	Porcine Relaxin (RLN) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 31.2-2000 pg/mL. Sensitivity: 8.36 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
<b>Expiration Date</b>	Please enquire.