

Product Datasheet

Rat Regenerating Islet Derived Protein 3 Alpha (REG3a) ELISA Kit (orb1663797)

Description	Rat Regenerating Islet Derived Protein 3 Alpha (REG3a) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 pg/mL. Sensitivity: 2.34 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.