

Product Datasheet

Mouse Regenerating Islet Derived Protein 3 Alpha (REG3a) ELISA Kit (orb1663798)

Description	Mouse Regenerating Islet Derived Protein 3 Alpha (REG3a) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.