

Product Datasheet

Porcine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit (orb1663822)

Description	Porcine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 1.15 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.