

Product Datasheet

Human RAB37, Member RAS Oncogene Family (RAB37) ELISA Kit (orb1663847)

Description	Human RAB37, Member RAS Oncogene Family (RAB37) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.095 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.