

Product Datasheet

Rat Prostaglandin Endoperoxide Synthase 2 (PTGS2) ELISA Kit (orb1663891)

Description	Rat Prostaglandin Endoperoxide Synthase 2 (PTGS2) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 25.7 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.