

## Product Datasheet

### Human presepsin ELISA Kit (orb1663959)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human presepsin ELISA Kit  |
| <b>Reactivity</b>          | Human  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detect Range: 31.2-2000 pg/mL.<br>Sensitivity: 11.64 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.<br>Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.  |