

Product Datasheet

Mouse Peroxiredoxin 4 (PRDX4) ELISA Kit (orb1663964)

Description Mouse Peroxiredoxin 4 (PRDX4) ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL. Sensitivity: 0.108 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.

Expiration Date Please enquire.