

Product Datasheet

Mouse Proopiomelanocortin (POMC) ELISA Kit (orb1664018)

| | |
|----------------------------|---|
| Description | Mouse Proopiomelanocortin (POMC) ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.156-10 ng/mL. Sensitivity: 0.045 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids. |
| Expiration Date | Please enquire. |