

Product Datasheet

Human Phospholipase A2, Group V (PLA2G5) ELISA Kit (orb1664065)

Description	Human Phospholipase A2, Group V (PLA2G5) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.105 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.