

Product Datasheet

Mouse Pyruvate Kinase, Muscle (PKM2) ELISA Kit (orb1664084)

Description	Mouse Pyruvate Kinase, Muscle (PKM2) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.144 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.