

Product Datasheet

General Prostaglandin F2 Alpha (PGF2a) ELISA Kit (orb1664150)

Description	General Prostaglandin F2 Alpha (PGF2a) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 11.67 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.