

Product Datasheet

Rat Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit (orb1664212)

Description	Rat Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 12.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.