

Product Datasheet

Mouse Programmed Cell Death Protein 1 Ligand 1 (PDCD1LG1) ELISA Kit (orb1664234)

Description	Mouse Programmed Cell Death Protein 1 Ligand 1 (PDCD1LG1) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.058 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.