

## Product Datasheet

# General Platelet Activating Factor (PAF) ELISA Kit (orb1664288)

<b>Description</b>	General Platelet Activating Factor (PAF) ELISA Kit
<b>Reactivity</b>	All
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.468-30 ng/mL. Sensitivity: 0.224 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
<b>Expiration Date</b>	Please enquire.