

Product Datasheet

Human Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit (orb1664302)

Description	Human Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.063 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.