

Product Datasheet

Porcine Osteopontin (OPN) ELISA Kit (orb1664329)

Description	Porcine Osteopontin (OPN) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5-4000 pg/mL. Sensitivity: 29.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.